

IIP45 antibody



Catalog Number: 111638

Product name

IIP45 antibody

Immunogen

[Human IIP45 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

IIP45 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of IIP45 recombinant protein (Accession Number: NM_021933). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111638(IIP45 Antibody) at dilution of 1:600