

# IIP45 antibody



Catalog Number: 111638

## Product name

IIP45 antibody

## Immunogen

[Human IIP45 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

IIP45 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue. Observed molecular weight by Western-blot: 45 kDa

## Preparation

This antibody was obtained by immunization of IIP45 recombinant protein (Accession Number: NM\_021933). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111638(IIP45 Antibody) at dilution of 1:600