

Rat Cadherin-8 / CDH8 (His Tag) recombinant protein



Catalog Number: 502638

General Information

Gene Name Synonym

Cadherin-8

Protein Construction

A DNA sequence encoding the rat CDH8 (NP_445845.2) (Met1-Met621) was expressed with a polyhistidine tag at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 30

Molecular Mass

The recombinant rat CDH8 comprises 603 amino

acids and predicts a molecular mass of 66.5 kDa. The apparent molecular mass of the recombinant protein is approximately 75-85 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat Cadherin-8 / CDH8 Protein (His Tag) SDS-PAGE