

EPHA5 Monoclonal Antibody



Catalog Number: 166638

Product name

EPHA5 Monoclonal Antibody

Specificity

Human

Antibody description

Monoclonal antibody to EPHA5

Preparation

Antigen: Recombinant protein of human EPHA5

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:10 - 1:50

Validations



Western blot - EPHA5 Monoclonal Antibody

Western blot analysis of extracts of U-251MG cells, using EPHA5 antibody. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.