

Anti-Dppa2 antibody



Catalog Number: 176638

Product name

Mouse IgG1

Anti-Dppa2 antibody

Applications

Specificity

WB

Mouse

Dilutions

Antibody description

WB: 1:500-1:2,000

Mouse monoclonal antibody to Dppa2

Preparation

This antigen of this antibody was peptide

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.

Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Validations




Fig1:: Western blot analysis of DPPA2 on F9 cell lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody (1/500) was used in 5% BSA at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1:5,000 dilution was used for 1 hour at room temperature.