

P. falciparum pfHPRT (His Tag) recombinant protein



Catalog Number: 502338

General Information

Protein Construction

A DNA sequence encoding the pfHPRT (XP_001349989.1) (Met1-Leu231) was expressed with a polyhistidine tag at the N-terminus.

Organism

P. falciparum

Expression Host

E. coli

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant pfHPRT consists 249 amino acids and predicts a molecular mass of 28.6 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Plasmodium falciparum (strain 3D7) HPRT Protein (His Tag)