

OLFM1 antibody

Catalog Number: 113338

Product name

OLFM1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

OLFM1 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue.

Preparation

This antibody was obtained by immunization of OLFM1 recombinant protein (Accession Number: NM_058199). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:113338(AMY antibody) at dilution of 1:50 (under 10x lens)