

# FNDC5 antibody

Catalog Number: 110738

## Product name

FNDC5 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

FNDC5 Rabbit Polyclonal antibody. Positive IHC detected in human skeletal muscle tissue. Positive IP detected in mouse liver tissue. Positive WB detected in mouse placenta tissue, mouse brain tissue, mouse liver tissue. Observed molecular weight by Western-blot: 25kDa - 30 kDa

## Preparation

This antibody was obtained by immunization of FNDC5 recombinant protein (Accession Number: NM\_001171941). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

IP: 1:200-1:1000

IHC: 1:20-1:200

## Validations



mouse liver, mouse skeletal muscle and rat liver were subjected to SDS PAGE followed by western blot with Catalog No:110738(FNDC5 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours



Immunohistochemistry of paraffin-embedded human skeletal muscle slide using Catalog No:110738(FNDC5 Antibody) at dilution of 1:50