

# SOX14 Polyclonal Antibody



Catalog Number: 161738

## Product name

SOX14 Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to SOX14

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 91-240 of human SOX14 (NP\_004180.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

## Validations



Immunohistochemistry - SOX14 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat liver using SOX14 antibody at dilution of 1:100



(40x lens).

Immunofluorescence - SOX14 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using SOX14 antibody .