

Anti-STAT1 Rabbit antibody



Catalog Number: 166738

Product name

Anti-STAT1 Rabbit antibody

Specificity

Human, Mouse

Antibody description

STAT1 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human STAT1

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IP, FC

Dilutions

WB: 1/1000-1/2000

IP: 1/50

FC: 1/50

Validations

