

Anti-StAR Rabbit antibody



Catalog Number: 168738

Product name

Anti-StAR Rabbit antibody

Specificity

Human

Antibody description

StAR Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human StAR

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

1:2000

1:50

Validations



Western blot detection of StAR in K562,3T3,HeLa cell lysates using StAR Rabbit mAb(1:1000 diluted).Predicted band size:32kDa.Observed band size:32kDa.