

# NCS1 antibody

Catalog Number: 113048

## Product name

NCS1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NCS1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, A431 cells, HepG2 cells, human brain tissue, mouse cerebellum tissue, PC-3 cells, rat brain tissue, SH-SY5Y cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 22 kDa

## Preparation

This antibody was obtained by immunization of NCS1 recombinant protein (Accession Number: NM\_014286). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:500-1:5000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113048(FREQ Antibody) at dilution of 1:2000



IP Result of anti-FREQ (IP:Catalog No:113048, 4ug; Detection:Catalog No:113048 1:1000) with mouse brain tissue lysate 4000ug.