

Anti-KCNC2 antibody

Catalog Number: 176048

Product name

Anti-KCNC2 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to KCNC2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human potassium channel kv3.2b. aa range:589-638

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.: Western Blot analysis of various cells using KV3.2 Polyclonal Antibody