

Anti-CSNK2A1 Rabbit antibody



Catalog Number: 168048

Product name

Anti-CSNK2A1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

CSNK2A1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human CSNK2A1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC


Dilutions

WB: 1/1000

IHC: 1/20

FC: 1/20

Validations



Western blot detection of CSNK2A1 in HeLa,CHO-K1,C6,mouse Brain,Ramos using CSNK2A1 Rabbit mAb(1:1000 diluted)