

BAF250 Polyclonal Antibody



Catalog Number: 165448

Product name

BAF250 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to BAF250

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1300-1400 of mouse BAF250 (NP_001074288.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:100

Validations



Immunofluorescence - BAF250 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using ARID1A antibody .