

MRPS21 antibody



Catalog Number: 112848

Product name

MRPS21 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MRPS21 Rabbit Polyclonal antibody. Positive IHC detected in human kidney.

Preparation

This antibody was obtained by immunization of MRPS21 recombinant protein (Accession Number: NM_018997). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human kidney using Catalog No:112848(MRPS21 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:112848(MRPS21 antibody) at dilution of 1:100 (under 40x lens)