

PKC delta antibody

Catalog Number: 113848

Product name

PKC delta antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PKC delta Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 70-80 kDa

Preparation

This antibody was obtained by immunization of PKC delta recombinant protein (Accession Number: XM_047448564). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:1000-1:10000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113848 (PRKCD antibody) at dilution of 1:2000



IP Result of anti-PRKCD (IP: Catalog No:113848, 4ug; Detection: Catalog No:113848 1:2000) with mouse brain tissue lysate 4000ug.