

# TAOK2 antibody

Catalog Number: 115848

## Product name

TAOK2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TAOK2 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in mouse brain tissue, HEK-293 cells, L02 cells. Observed molecular weight by Western-blot: 138-150 kDa

## Preparation

This antibody was obtained by immunization of TAOK2 recombinant protein (Accession Number: NM\_016151). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115848(TAOK2 antibody) at dilution of 1:300



IP Result of anti-TAOK2 (IP:Catalog No:115848, 4ug; Detection:Catalog No:115848 1:500) with HEK-293 cells lysate 3000ug.