

Anti-SFPQ Rabbit antibody



Catalog Number: 166848

Product name

Anti-SFPQ Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

SFPQ Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human SFPQ

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/2000-1/10000

IHC: 1/50

ICC/IF: 1/20

Validations



Western blot detection of SFPQ in K562,Rat Brain,C6,3T3,Hela cell lysates using SFPQ Rabbit pAb(1:1000 diluted).Predicted band size:76kDa.Observed band size:100kDa.