

# Iladatuzumab vedotin reference antibody



Catalog Number: 176848

## Product name

Iladatuzumab vedotin reference antibody

## Specificity

Human CD79b

## Antibody description

Anti-CD79b Antibody (Iladatuzumab vedotin)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations




## SDS-PAGE

Anti-CD79b Antibody (iladatuzumab vedotin) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-CD79b Antibody (Iladatuzumab vedotin) is more than 95%, determined by SEC-HPLC.



## HIC-HPLC

The drug-to-antibody ratio (DAR) of the Anti-CD79b Antibody (iladatuzumab vedotin) is 3.8.

# Iladatuzumab vedotin reference antibody



Catalog Number: 176848

## ELISA

Immobilized human CD79b His at 2 ug/mL can bind Anti-CD79b Antibody (iladatuzumab vedotin), EC50 is 0.01675 ug/mL.

