

# Anti-HMGB1 Rabbit antibody



Catalog Number: 168848

## Product name

Anti-HMGB1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

HMGB1 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human HMGB1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions


1:2000

1:50

1:50

1:20

## Validations



Western blot detection of HMGB1 in K562, Rat Brain, C6, 3T3, HeLa cell lysates using HMGB1 Rabbit mAb(1:1000 diluted). Predicted band size: 25kDa. Observed band size: 25kDa.