

# APAF1 Polyclonal Antibody



Catalog Number: 160148

## Product name

APAF1 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to APAF1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human APAF1 (NP\_863651.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - APAF1 Polyclonal Antibody

Western blot analysis of extracts of 293 cells, using APAF1 antibody .