

Anti-OR10H3 antibody



Catalog Number: 176148

Product name

Polyclonal

Anti-OR10H3 antibody

Ig Type

Specificity

Rabbit IgG

Human

Applications

Antibody description

WB

Rabbit polyclonal antibody to OR10H3

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or10h3/10h4. aa range:234-283

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Validations




Fig1:: Western Blot analysis of HELA using Olfactory receptor 10H3/4 Polyclonal Antibody. Antibody was diluted at 1:1000