

HLA-DOA Polyclonal Antibody



Catalog Number: 162548

Product name

HLA-DOA Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to HLA-DOA

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 26-217 of human HLA-DOA (NP_002110.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:200 - 1:2000

Validations



Western blot - HLA-DOA Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HLA-DOA antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90S.