

p35 antibody

Catalog Number: 113548

Product name

p35 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

p35 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of p35 recombinant protein (Accession Number: XM_011525438). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113548(CDK5R1 antibody) at dilution of 1:600



IP Result of anti-CDK5R1 (IP:Catalog No:113548, 4ug; Detection:Catalog No:113548 1:500) with mouse brain tissue lysate 3440ug.