

# SPATA13 antibody



Catalog Number: 115548

## Product name

SPATA13 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

SPATA13 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 75kd

## Preparation

This antibody was obtained by immunization of SPATA13 recombinant protein (Accession Number: NM\_153023). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:115548(SPATA13 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using SPATA13 antibody Catalog No:115548 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).