

Phospho-PERK-T982 pAb



Catalog Number: 166548

Product name

Phospho-PERK-T982 pAb

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to Phospho-PERK-T982

Preparation

Antigen: A phospho specific peptide corresponding to residues surrounding T980 of human PERK.

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-PERK-T982 pAb

Western blot analysis of extracts of C6 cells, using Phospho-PERK-T982 pAb at 1:1000 dilution. C6 cells were treated by ATP (5 mM) at 30u2103 for 1 hour. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA. Detection: ECL Basic Kit. Exposure time: 1s.