

Anti-Caspase-9 Rabbit antibody



Catalog Number: 167548

Product name

Anti-Caspase-9 Rabbit antibody

Specificity

Mouse, Rat

Antibody description

Caspase-9 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of mouse Caspase-9

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions


WB: 1/1000

IP: 1/20

Validations



Western blot detection of Caspase-9 in 3T3 cell lysates using Caspase-9 Rabbit pAb(1:1000 diluted). Predicted band size:50kDa. Observed band



size:50kDa.

Western blot detection of Caspase-9 in C6 cell lysates using Caspase-9 Rabbit pAb(1:1000 diluted). Predicted band size:50kDa. Observed band size:50kDa.