

Human PPIA/CYPA recombinant protein



Catalog Number: 518548

General Information

Gene Name Synonym

Peptidyl-prolyl cis-trans isomerase A; PPIase A; Cyclophilin A; Cyclosporin A-binding protein; Rotamase A; SP18; PPIA; CYPA

Protein Construction

A DNA sequence encoding the Human PPIA/CYPA Met1-Glu165 is expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of PBS, 10% glycerol, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PPIA/CYPA recombinant protein