

SERPINA12 Polyclonal Antibody



Catalog Number: 161948

Product name

SERPINA12 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to SERPINA12

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-320 of human SERPINA12 (NP_776249.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC, IF

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200

Validations

Immunofluorescence - SERPINA12 Polyclonal Antibody

Immunofluorescence analysis of PC-12 cells using SERPINA12 antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.