

Anti-DDX55 antibody

Catalog Number: 175948

Product name

Anti-DDX55 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to DDX55

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human ddx55. aa range:91-140

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of various cells using DDX55 Polyclonal Antibody diluted at 1: 500