

ZNF140 antibody



Catalog Number: 116948

Product name

ZNF140 antibody

Immunogen

[Human ZNF140 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ZNF140 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in SKOV-3 cells, Jurkat cells. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of ZNF140 recombinant protein (Accession Number: NM_003440). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



SKOV-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116948(ZNF140 Antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells using Catalog No:116948(ZNF140 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG