

# CBX6 antibody

Catalog Number: 108948

## Product name

CBX6 antibody

## Specificity

Human; other species not tested.

## Antibody description

CBX6 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Positive IHC detected in human skin tissue. Observed molecular weight by Western-blot: 55-60 kDa

## Preparation

This antibody was obtained by immunization of CBX6 recombinant protein (Accession Number: NM\_014292). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:108948(CBX6 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human skin slide using Catalog No:108948(CBX6 Antibody) at dilution of 1:50