

Endoglin/CD105 antibody

Catalog Number: 110248

Product name

Endoglin/CD105 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Endoglin/CD105 Rabbit Polyclonal antibody. Positive IHC detected in human placenta tissue. Positive IP detected in mouse lung tissue. Positive WB detected in mouse lung tissue. Observed molecular weight by Western-blot: 170-180 kDa

Preparation

This antibody was obtained by immunization of Endoglin/CD105 recombinant protein (Accession Number: NM_001114753). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IP, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:1000


Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:110248(ENG antibody) at dilution of 1:300



IP Result of anti-ENG (IP:Catalog No:110248, 3ug; Detection:Catalog No:110248 1:300) with mouse lung tissue lysate 4000ug.




Immunohistochemical of paraffin-embedded human placenta using Catalog No:110248(ENG

Endoglin/CD105 antibody

Catalog Number: 110248

antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human placenta using Catalog No:110248(ENG antibody) at dilution of 1:50 (under 40x lens)