

Anti-MID1IP1 antibody



Catalog Number: 141248

Product name

Anti-MID1IP1 antibody

Immunogen

[Human MID1IP1 \(His Tag\) recombinant protein](#)

Specificity

Human MID1IP1

Antibody description

Rabbit Monoclonal to Human MID1IP1

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human MID1IP1 .

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 5-20 µg/ml

Validations



Lanes	A
Items	
Sample (whole cell lysate)	HepG2
Sample Volume (µg/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	5-20 µg/ml
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.
Developed using Odyssey imaging system.	
Explanation	Predicted band size : 20 kDa Observed band size : 23 kDa

Human MID1IP1 Western blot (WB) 14852