

PAEP Polyclonal Antibody



Catalog Number: 162248

Product name

PAEP Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PAEP

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 19-180 of human PAEP (NP_002562.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:200

Validations



Immunohistochemistry - PAEP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using PAEP antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - PAEP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse pancreas using PAEP antibody at dilution of 1:100 (40x lens).