

Anti-Bok Rabbit antibody



Catalog Number: 167248

Product name

Anti-Bok Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Bok Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Bok

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, FC

Dilutions

WB: 1/1000

FC: 1/20

Validations



Western blot detection of Bok in Rat Brain,C6,HeLa cell lysates using Bok Rabbit pAb(1:1000 diluted).Predicted band size:23kDa.Observed band size:23kDa.