

CXorf15 antibody



Catalog Number: 109648

Product name

CXorf15 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CXorf15 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HEK-293 cells, Jurkat cells, MCF-7 cells. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 66 kDa

Preparation

This antibody was obtained by immunization of CXorf15 recombinant protein (Accession Number: NM_018360). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:109648(CXorf15 Antibody) at dilution of



Immunofluorescent analysis of HEK-293 cells using Catalog No:109648(CXorf15 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)