

Anti-THAP4 antibody

Catalog Number: 176348

Product name

Anti-THAP4 antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

Rabbit polyclonal antibody to THAP4

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human thap4. aa range:341-390

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.: Western Blot analysis of various cells using THAP4 Polyclonal Antibody