

# Trappc9,NIBP antibody



Catalog Number: 116348

## Product name

Polyclonal

Trappc9,NIBP antibody

## Ig Type

## Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

Trappc9,NIBP Rabbit Polyclonal antibody. Positive WB detected in Neuro-2a cells, mouse brain tissue, mouse skeletal muscle tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 128kd,139kd

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

## Validations

This antibody was obtained by immunization of Trappc9,NIBP recombinant protein (Accession Number: NM\_180662). Purification method: Antigen affinity purified.



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Neuro-2a cells were subjected to SDS PAGE followed by western blot with Catalog No:116348(Trappc9,NIBP antibody) at dilution of 1:300