

INHBA-Specific antibody



Catalog Number: 107348

Product name

INHBA-Specific antibody

Specificity

Human, Pig; other species not tested.

Antibody description

INHBA-Specific Mouse Monoclonal antibody. Positive WB detected in human testis tissue, pig liver tissue. Positive IHC detected in human lung cancer tissue, human oesophagus cancer tissue, human stomach cancer tissue. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002192). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:500-1:5000


IHC: 1:50-1:500

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107348(INHBA-Specific Antibody) at dilution of 1:1000

Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog No:107348(INHBA-Specific Antibody) at dilution of 1:200 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog No:107348(INHBA-Specific Antibody) at dilution of 1:150 (under 40x lens). Heat mediated antigen retrieved with Tris-EDTA buffer(pH9).