

# CLEC14A antibody



Catalog Number: 109348

## Product name

CLEC14A antibody

## Immunogen

[Human CLEC14A Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CLEC14A Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 52 kDa

## Preparation

This antibody was obtained by immunization of CLEC14A recombinant protein (Accession Number: NM\_175060). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109348(CLEC14A antibody) at dilution of 1:1000