

Human Moesin (aa 1-346, His Tag) recombinant protein



Catalog Number: 500748

General Information

Gene Name Synonym

Membrane-organizing extension spike protein

Protein Construction

A DNA sequence encoding the human MSN (P26038) (Met 1-Glu 346) was expressed, with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 80 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human MSN consists of 361 amino acids and has a calculated molecular mass

of 42.8 kDa. It migrates as an approximately 45 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 0.5M NaCl, pH 8.0

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human MSN / Moesin Protein (aa 1-346, His Tag)
SDS-PAGE