

# IL-6 antibody

Catalog Number: 111748

## Product name

IL-6 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

IL-6 Rabbit Polyclonal antibody. Positive WB detected in Recombinant protein, Transfected HEK-293 cells. Positive IHC detected in human tonsillitis tissue, human spleen tissue. Observed molecular weight by Western-blot: 30 kDa

## Preparation

This antibody was obtained by immunization of IL-6 recombinant protein (Accession Number: XM\_011515390). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:111748(IL6 Antibody) at dilution of 1:3000



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111748(IL6 Antibody) at dilution of 1:1000




Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111748(IL6 Antibody) at dilution of 1:200 (under 10x lens)

## IL-6 antibody

Catalog Number: 111748

---



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111748(IL6 Antibody) at dilution of 1:200 (under 40x lens)