

SETD1B Polyclonal Antibody



Catalog Number: 163748

Product name

SETD1B Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to SETD1B

Preparation

Antigen: Recombinant protein of human SETD1B

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF, IP

Dilutions


WB 1:200 - 1:2000

IHC 1:50- 1:200

IF 1:50- 1:200

IP 1:20- 1:50

Validations



Western blot - SETD1B Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SETD1B antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.