

DCD antibody

Catalog Number: 109748

Product name

DCD antibody

Specificity

Human; other species not tested.

Antibody description

DCD Rabbit Polyclonal antibody. Positive IHC detected in human skin tissue, human breast cancer tissue.

Preparation

This antibody was obtained by immunization of DCD recombinant protein (Accession Number: NM_053283). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:109748(DCD Antibody) at dilution of 1:800 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:109748(DCD Antibody) at dilution of 1:800 (under 40x lens)