

# iNOS Polyclonal Antibody



Catalog Number: 160058

## Product name

iNOS Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to iNOS

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 500-730 of human iNOS (NP\_000616.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



## Western blot - iNOS Polyclonal Antibody

Western blot analysis of extracts of Raw 264.7 cells, using iNOS antibody at 1:1000 dilution. Raw264.7 cells were treated by LPS (1 µg/ml) for 8 hours. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.