

# GPM6A Polyclonal Antibody



Catalog Number: 166058

## Product name

GPM6A Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to GPM6A

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GPM6A (NP\_005268.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

## Validations




Western blot - GPM6A Polyclonal Antibody

Western blot analysis of extracts of U-87MG cells, using GPM6A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 30s.



Immunohistochemistry - GPM6A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse stomach using GPM6A antibody at dilution



of 1:100 (40x lens).

# GPM6A Polyclonal Antibody



Catalog Number: 166058

---

Immunohistochemistry - GPM6A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using GPM6A antibody at dilution of 1:100 (40x lens).