

THTPA antibody

Catalog Number: 116058

Product name

THTPA antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

THTPA Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in K-562 cells, A431 cells, HeLa cells, human brain tissue, mouse testis tissue, mouse uterus tissue, PC-3 cells. Positive IF detected in HepG2 cells. Positive IHC detected in human stomach cancer tissue. Observed molecular weight by Western-blot: 26-35kd

Preparation

This antibody was obtained by immunization of THTPA recombinant protein (Accession Number: NM_024328). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:10-1:100

Validations




K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:116058(THTPA antibody) at dilution of 1:800




Immunohistochemical of paraffin-embedded human stomach cancer using Catalog No:116058(THTPA antibody) at dilution of 1:50 (under 40x lens)

THTPA antibody

Catalog Number: 116058



Immunofluorescent analysis of HepG2 cells, using THTPA antibody Catalog No:116058 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-THTPA (IP:Catalog No:116058, 3ug; Detection:Catalog No:116058 1:800) with K-562 cells lysate 1720ug.