

# Anti-PDLIM7 Rabbit antibody



Catalog Number: 167058

## Product name

Rabbit IgG

Anti-PDLIM7 Rabbit antibody

## Applications

## Specificity

WB, IHC-P, FC

Human

## Dilutions

## Antibody description

WB: 1/1000

PDLIM7 Rabbit polyclonal antibody

IHC: 1/20

## Preparation


FC: 1/20

Antigen: Recombinant protein of human PDLIM7

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



## Storage

Store at -20°C. Stable for 12 months from date of receipt.

Western blot detection of PDLIM7 in HeLa cell lysates using PDLIM7 Rabbit pAb(1:1000 diluted).Predicted band size:50kDa.Observed band size:50kDa.

## Clonality

Polyclonal

## Ig Type